Surface Plasmon Resonance Sensors: Optimization of Diffraction Grating and Prism Couplers

W. Raja¹, A. Alabastri¹, S. Tuccio¹, R. Proietti Zaccaria¹

Abstract

Surface plasmon resonance (SPR) sensors proved themselves as a promising device for many kinds of applications such as optical biosensing [1], binding constant determination [2] or nanofilm thickness measurements [3]. Here we simulate using COMSOL Multiphysics® the light-polaritons coupling for the two most commonly used SPR setups: Attenuated total reflection (Kretschmann configuration) and diffraction gratings as shown in Figure 1. Furthermore, we illustrate a comprehensive designing rule to determine the parameters for optimizing the generation of surface plasmon polaritons (SPP) for the both setups i.e. grating couplers and Kretschmann configuration. Moreover, a computational analysis are demonstrated for three detection methods: Resonant angle (Figure 2), Resonant wavelength (Figure 3) and Thickness of metal (Figure 4), in SPR sensors for both the SPR setups.

Reference

- [1] Jiří Homola, Sinclair S. Yee and Günter Gauglitz, Sensors and Actuators B: Chemical 54, 3 (1999).
- [2] Matthias C. Jecklin, Stefan Schauer, Christoph E. Dumelinand Renato Zenobi, J. Mol. Recognit.22, 319 (2009).
- [3] Guillaume Legay, Eric Finot, Rita Meunier-Prest, Mustapha Cherkaoui-Malki, Norbert Latruffe, Alain Dereux, Biosensors and Bioelectronics 21, 627 (2005).
- [4] David B. Hand, J. Biol. Chem. 108:703-707 (1935).
- [5] S. Lousinian, S. Logothetidis Microelectronic Engineering 84, 479–485 (2007).

¹Department of Nanostructures, Istituto Italiano di Tecnologia, Genova, Italy

Figures used in the abstract




Figure 1: (a) Surface Plasmon Polariton (SPP) generation for a transverse magnetic mode (TM). XY: plane of incidence; Z: polarization direction; X: direction of propagation of the generated SPP. The Hz component is shown. (b) Surface Plasmon Resonance through a grating coupler with 23 grating element with a period of 580nm. The Hz component is shown (c). Surface plasmon Resonance demonstration through single element with Floquet periodic boundary condition with a period of 580nm. The norm of the magnetic field H is shown. R.I. stands for Refractive Index.




Figure 2: (a) Reflectivity as a function of the angle of incidence calculated at 632nm (Helium Neon Laser). Various proteins (refractive indexes) taken from human and animal blood samples were considered [4] in Kretschmann configuration. (b) Reflectivity and absorption as a function of the angle of incidence when fibronogen (FIB) was considered [5] in grating couplers.




Figure 3: (a) Reflectivity as a function of wavelength calculated for various proteins taken from human and animal blood sample [4] in Kretschmann configuration. (b) Reflectivity as a function of wavelength human serum Albumin taken from human blood sample [5] in grating couplers.




Figure 4: Reflectivity as a function of the thickness of the metal (d) with respect to the incidence angle calculated for the fibronogen protein (FIB) taken from human blood sample [5] in in Kretschmann configuration.